

[最近のトピックス]

Measurement of IP₃ by the Competitive Fluorescent Ligand Assay (CFLA) method

Azmerce JAHAN, Akihiro NEZU, and Akihiko TANIMURA

Division of Pharmacology, School of Dentistry, Health Sciences University of Hokkaido

Key words : Inositol 1,4,5-trisphosphate, Fluorescence Resonance Energy Transfer, Fluorescent sensor

Inositol 1,4,5-trisphosphate (IP₃) is an important intracellular messenger in Ca²⁺ signaling and is involved in numerous responses to hormones, neurotransmitters, and growth factors. The IP₃ releases Ca²⁺ from intracellular stores through the IP₃ receptor (IP₃R), and this makes IP₃ measurement advantageous. There are several methods to measure IP₃, including (1) AlphaScreen technology (Taouji et al., 2009), and (2) fluorescence polarization for detection of the binding of IP₃ binding proteins and fluorescent analogues of IP₃ (Rossi & Taylor, 2011). These methods, however, required large amounts of IP₃ binding protein and accurate concentrations of functional binding sites. Presently, a conceptually new method for measuring IP₃, the Competitive Fluorescent Ligand Assay for IP₃ (CFLA-IP₃) has become available.

The CFLA-IP₃ uses Fluorescence Resonance Energy Transfer (FRET) which occurs between two different fluorescent molecules, CFP-labeled ligand binding domains of IP₃R (CFP-LBD), and fluorescent ligand (FL). The binding of FL to CFP-LBD causes FRET, and IP₃ competes with the FL in the binding to the LBD, and thereby decreases the FRET signal (Fig. 1).

Oura et al., developed a new high affinity fluorescent ligand for IP₃R, fluorescent adenophostin A by introducing the fluorescein unit at the 5'-position of ADA and its low affinity analogue (F-LL). It has been shown that the binding of these FL and CFP-LBD decreased the CFP signal and increased the FL signal due to FRET. In addition, F-LL-induced changes in the fluorescence ratio (CFP/FL) were reduced by the addition of IP₃ in a concentration-dependent manner. This method was further applied to measure the IP₃ concentration in cytosolic fractions of COS-7 cells with and without stimulation by ATP, and to examine the potency of

IP₃R ligands.

These results show that the CFLA method is reliable, and the principles of this method are simply applicable to any receptor for which a fluorescent ligand is available.

References

- Oura T, Murata K, Morita T, Nezu A, Arisawa M, Shuto S & Tanimura A. Highly sensitive measurement of Inositol 1,4,5-trisphosphate by using a new fluorescent ligand and ligand binding domain combination. *Chembiochem* 17 : 1509–1512, 2016.
- Rossi AM & Taylor CW. Analysis of protein–ligand interactions by fluorescence polarization. *Nat Protoc* 6 : 365–387, 2011.
- Taouji S, Dahan S, Bossé R & Chevet E. Current Screens Based on the AlphaScreen Technology for Deciphering Cell Signalling Pathways. *Curr Genomics* 10 : 93–101, 2009.




Figure 1. Competitive Fluorescent Ligand Assay for IP₃ (CFLA-IP₃). FL : Fluorescent ligand, CFP : Cyan fluorescent protein, LBD : Ligand binding domain of IP₃R.